

Product datasheet for **TA398548**

MYOD1 Rabbit Monoclonal Antibody [Clone ID: RM369]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM369
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to residues near N-terminus of human MyoD1 (Myoblast determination protein 1)
Specificity:	This antibody reacts to human, mouse, and rat MyoD1 (Myoblast determination protein 1)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	myogenic differentiation 1
Database Link:	Entrez Gene 4654 Human P15172
Synonyms:	bHLHc1; Myf-3; MYF3; MYOD; PUM



[View online »](#)