

Product datasheet for **TA398549**

SFTP3 Rabbit Monoclonal Antibody [Clone ID: RM370]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM370
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:500-1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to residues near C-terminus of human SP-B (Surfactant protein B)
Specificity:	This antibody reacts to human SP-B (Surfactant protein B) and proSP-B
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	surfactant protein B
Database Link:	Entrez Gene 6439 Human P07988
Synonyms:	OTTHUMP00000203253; PSP-B; SFTB3; SFTP3; SMDP1; SP-B



[View online »](#)