

Product datasheet for **TA398553**

c Fos (FOS) Rabbit Monoclonal Antibody [Clone ID: RM374]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM374
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:1000-1:5000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to N-terminus of human Proto-oncogene c-Fos
Specificity:	This antibody reacts to human Proto-oncogene c-Fos. It may also react to mouse and rat c-Fos, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	Fos proto-oncogene, AP-1 transcription factor subunit
Database Link:	Entrez Gene 2353 Human P01100
Synonyms:	AP-1; C-FOS; G0S7



[View online »](#)