

Product datasheet for **TA398554**

MSH2 Rabbit Monoclonal Antibody [Clone ID: RM375]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM375
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:200-1:1000 dilution WB (WB): 1:400-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A mammalian cells expressed recombinant protein fragment within 500 amino acids of the N-terminus of human Msh2
Specificity:	This antibody reacts to human DNA mismatch repair protein Msh2
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	mutS homolog 2
Database Link:	Entrez Gene 4436 Human P43246
Synonyms:	COCA1; FCC1; HNPCC; HNPCC1; LCFS2



[View online »](#)