

Product datasheet for **TA398555**

MSH6 Rabbit Monoclonal Antibody [Clone ID: RM376]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM376
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:500 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A mammalian cells expressed recombinant protein fragment of human Msh6
Specificity:	This antibody reacts to human DNA mismatch repair protein Msh6
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	mutS homolog 6
Database Link:	Entrez Gene 2956 Human P52701
Synonyms:	GTBP; GTMBP; HNPCC5; HSAP; p160



[View online »](#)