

Product datasheet for **TA398556**

Arginase 1 (ARG1) Rabbit Monoclonal Antibody [Clone ID: RM377]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM377
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:1000-1:2000 dilution WB (WB): 1:5000-1:10000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human Arginase-1 (ARG1)
Specificity:	This antibody reacts to human Arginase-1 (ARG1)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	arginase 1
Database Link:	Entrez Gene 383 Human P05089
Synonyms:	ARG1



[View online »](#)