

## Product datasheet for **TA398564**

### **CD61 (ITGB3) Rabbit Monoclonal Antibody [Clone ID: RM382]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	RM382
<b>Applications:</b>	IHC, WB
<b>Recommended Dilution:</b>	Immunohistochemistry (IHC): 1:250-1:500 dilution WB (WB): 1:5000-1:25000 dilution
<b>Reactivity:</b>	Human
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	A peptide corresponding to the internal region of human Integrin $\beta$ 3.
<b>Specificity:</b>	This antibody reacts to Human Integrin beta3A / CD61
<b>Formulation:</b>	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
<b>Concentration:</b>	N/A - lot specific
<b>Purification:</b>	Protein A affinity purified from an animal origin-free culture supernatant
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Gene Name:</b>	integrin subunit beta 3
<b>Database Link:</b>	<a href="#">Entrez Gene 3690 Human P05106</a>
<b>Synonyms:</b>	CD61; GP3A; GPIIIa



[View online »](#)