

## Product datasheet for **TA398567**

### Chromogranin A (CHGA) Rabbit Monoclonal Antibody [Clone ID: RM385]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM385
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:5000-1:10000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human Chromogranin A
Specificity:	This antibody reacts to Human Chromogranin A
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	chromogranin A
Database Link:	<a href="#">Entrez Gene 1113 Human P10645</a>
Synonyms:	CgA; SP-I



[View online »](#)