

Product datasheet for **TA398568**

Cytokeratin 10 (KRT10) Rabbit Monoclonal Antibody [Clone ID: RM386]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM386
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:200-1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human CK-10
Specificity:	This antibody reacts to Human CK-10
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	keratin 10
Database Link:	Entrez Gene 3858 Human P13645
Synonyms:	CK-10; CK10; Cytokeratin-10; K10; Keratin-10; KPP



[View online »](#)