

Product datasheet for **TA398571**

p63 (TP63) Rabbit Monoclonal Antibody [Clone ID: RM392]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM392
Applications:	IHC
Recommended Dilution:	Immunohistochemistry (IHC): 1:50-1:200 dilutio
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to Human p40 (Δ Np63)
Specificity:	This antibody reacts to human p40 (Δ Np63)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	tumor protein p63
Database Link:	Entrez Gene 8626 Human Q9H3D4
Synonyms:	AIS; B(p51A); B(p51B); CUSP; EEC3; KET; LMS; NBP; OFC8; OTTHUMP00000209738; OTTHUMP00000209740; p40; p51; p53CP; p63; p73H; p73L; RHS; SHFM4; TP53CP; TP53L; TP73L



[View online »](#)