

Product datasheet for **TA398572**

CD1 (CD1A) Rabbit Monoclonal Antibody [Clone ID: RM393]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM393
Applications:	IHC
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:200 dilutio
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to Human CD1a
Specificity:	This antibody reacts to human CD1a
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	CD1a molecule
Database Link:	Entrez Gene 909 Human P06126
Synonyms:	CD1; FCB6; HTA1; R4; T6



[View online »](#)