

Product datasheet for **TA398573**

Topoisomerase II alpha (TOP2A) Rabbit Monoclonal Antibody [Clone ID: RM394]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM394
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:2000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to Human TOP2a
Specificity:	This antibody reacts to human TOP2a
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	topoisomerase (DNA) II alpha
Database Link:	Entrez Gene 7153 Human P11388
Synonyms:	OTTHUMP00000207108; TOP2; TP2A



[View online »](#)