

Product datasheet for **TA398575**

Caldesmon (CALD1) Rabbit Monoclonal Antibody [Clone ID: RM396]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM396
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:200 dilution WB (WB): 1:100-1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to Human CALD1
Specificity:	This antibody reacts to human CALD1
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	caldesmon 1
Database Link:	Entrez Gene 800 Human Q05682
Synonyms:	CAD; CDM; H-CAD; HCAD; L-CAD; LCAD; MGC21352; NAG22; OTTHUMP00000207949; OTTHUMP00000207954



[View online »](#)