

## Product datasheet for **TA398578**

### VEGFA Rabbit Monoclonal Antibody [Clone ID: RM391]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM391
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:200 dilution WB (WB): 1:100-1:200 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to Human VEGFA
Specificity:	This antibody reacts to human VEGFA
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	vascular endothelial growth factor A
Database Link:	<a href="#">Entrez Gene 7422 Human P15692</a>
Synonyms:	MGC70609; MVCD1; VEGF; VEGF-A; VPF



[View online »](#)