

Product datasheet for **TA398579**

Rabbit Monoclonal Antibody [Clone ID: RM400]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM400
Applications:	IHC
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:200 dilutio
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Proteins purified from Helicobacter pylori
Specificity:	This antibody reacts to Helicobacter pylori
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Database Link:	=



[View online »](#)