

Product datasheet for **TA398580**

FOXP1 Rabbit Monoclonal Antibody [Clone ID: RM402]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM402
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:200 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to residues near the C-terminus of human Foxp1
Specificity:	This antibody reacts to human to Foxp1
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	forkhead box P1
Database Link:	Entrez Gene 27086 Human Q9H334
Synonyms:	12CC4; FLJ23741; hFKH1B; HSPC215; MGC12942; MGC88572; MGC99551; QRF1



[View online »](#)