

Product datasheet for **TA398585**

Myeloperoxidase (MPO) Rabbit Monoclonal Antibody [Clone ID: RM407]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM407
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:200 dilution WB (WB): 1:200-1:500 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the internal region of light chain on MPO. This antibody may also react to mouse or rat MPO, as predicted by immunogen homology.
Specificity:	This antibody reacts to human to MPO
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	myeloperoxidase
Database Link:	Entrez Gene 4353 Human P05164
Synonyms:	myeloperoxidase



[View online »](#)