

Product datasheet for **TA398589**

Cytokeratin 7 (KRT7) Rabbit Monoclonal Antibody [Clone ID: RM416]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM416
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 -1:200 dilution WB: 1:1000 - 1:10000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human Cytokeratin-7
Specificity:	This antibody reacts to human CK-7
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	keratin 7
Database Link:	Entrez Gene 3855 Human P08729
Synonyms:	CK-7; CK7; Cytokeratin-7; K2C7; K7; Keratin-7; MGC3625; MGC129731; sarcolectin; SCL



[View online »](#)