

Product datasheet for **TA398591**

IL2 Receptor alpha (IL2RA) Rabbit Monoclonal Antibody [Clone ID: RM418]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM418
Applications:	IHC
Recommended Dilution:	IHC: 1:100 -1:200 dilutio
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human CD25
Specificity:	This antibody reacts to human CD25 (IL2R α)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	interleukin 2 receptor subunit alpha
Database Link:	Entrez Gene 3559 Human P01589
Synonyms:	CD25; IDDM10; IL-2-RA; IL2-RA; IL2R; p55; TCGFR



[View online »](#)