

Product datasheet for **TA398594**

Annexin A1 (ANXA1) Rabbit Monoclonal Antibody [Clone ID: RM424]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM424
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 -1:200 dilution WB: 1:1000 - 1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human Annexin A1
Specificity:	This antibody reacts to human Annexin A1
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	annexin A1
Database Link:	Entrez Gene 301 Human P04083
Synonyms:	Annexin-1; ANX1; Chromobindin-9; LPC1; OTTHUMP00000021475; p35



[View online »](#)