

## Product datasheet for **TA398596**

### Prostein (SLC45A3) Rabbit Monoclonal Antibody [Clone ID: RM426]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM426
Applications:	IHC
Recommended Dilution:	IHC: 1:100 -1:400 dilutio
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human Prostein
Specificity:	This antibody reacts to human Prostein
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	solute carrier family 45 member 3
Database Link:	<a href="#">Entrez Gene 85414 Human</a> <a href="#">Q96JT2</a>
Synonyms:	IPCA-2; IPCA-6; IPCA-8; IPCA6; PCANAP2; PCANAP6; PCANAP8; prostein; PRST



[View online »](#)