

Product datasheet for **TA398599**

Nkx3.1 (NKX3-1) Rabbit Monoclonal Antibody [Clone ID: RM430]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM430
Applications:	IHC
Recommended Dilution:	IHC: 1:100 -1:400 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human Nkx3.1
Specificity:	This antibody reacts to human Nkx3.1
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	NK3 homeobox 1
Database Link:	Entrez Gene 4824 Human Q99801



[View online »](#)