

Product datasheet for **TA398601**

Rabbit Monoclonal Antibody [Clone ID: RM432]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM432
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 – 1:500 dilution WB: 1:2000 - 1:5000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A recombinant human HSP70 protein
Specificity:	This antibody reacts to human HSP70
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Database Link:	P0DMV8 , P0DMV9



[View online »](#)