

## Product datasheet for **TA398603**

### **PAX8 Rabbit Monoclonal Antibody [Clone ID: RM436]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	RM436
Applications:	IHC
Recommended Dilution:	IHC: 1:100 -1:400 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of PAX8
Specificity:	This antibody reacts to PAX8
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	paired box 8
Database Link:	<a href="#">Entrez Gene 7849 Human Q06710</a>
Synonyms:	PAX8



[View online »](#)