

## Product datasheet for **TA398604**

### Clusterin (CLU) Rabbit Monoclonal Antibody [Clone ID: RM437]

#### Product data:

|                       |                                                                                               |
|-----------------------|-----------------------------------------------------------------------------------------------|
| Product Type:         | Primary Antibodies                                                                            |
| Clone Name:           | RM437                                                                                         |
| Applications:         | IHC, WB                                                                                       |
| Recommended Dilution: | IHC: 1:100 -1:400 dilution<br>WB: 1:1000-1:2000 dilution                                      |
| Reactivity:           | Human                                                                                         |
| Host:                 | Rabbit                                                                                        |
| Isotype:              | IgG                                                                                           |
| Clonality:            | Monoclonal                                                                                    |
| Immunogen:            | A peptide corresponding to the C-terminus of Clusterin                                        |
| Specificity:          | This antibody reacts to Clusterin                                                             |
| Formulation:          | 50% Glycerol/PBS with 1% BSA and 0.09% sodium azide                                           |
| Concentration:        | N/A - lot specific                                                                            |
| Purification:         | Protein A affinity purified from an animal origin-free culture supernatant                    |
| Conjugation:          | Unconjugated                                                                                  |
| Storage:              | Store at -20°C as received.                                                                   |
| Stability:            | Stable for 12 months from date of receipt.                                                    |
| Gene Name:            | clusterin                                                                                     |
| Database Link:        | <a href="#">Entrez Gene 1191 Human P10909</a>                                                 |
| Synonyms:             | AAG4; Apo-J; APOJ; CLI; clusterin; KUB1; MGC24903; NA1/NA2; SGP-2; SGP2; SP-40; TRPM-2; TRPM2 |



[View online »](#)