

## Product datasheet for **TA398605**

### **S6K1 (RPS6KB1) Rabbit Monoclonal Antibody [Clone ID: RM438]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	RM438
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 – 1:400 dilution WB: 1:1000 - 1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of p70S6K
Specificity:	This antibody reacts to p70S6K
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ribosomal protein S6 kinase B1
Database Link:	<a href="#">Entrez Gene 6198 Human P23443</a>
Synonyms:	p70(S6K)-alpha; p70-alpha; p70-S6K; P70S6K; PS6K; S6K; S6K1; STK14A



[View online »](#)