

Product datasheet for **TA398606**

HDAC3 Rabbit Monoclonal Antibody [Clone ID: RM439]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM439
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 – 1:200 dilution WB: 1:1000 - 1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of HDAC3
Specificity:	This antibody reacts to HDAC3
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	histone deacetylase 3
Database Link:	Entrez Gene 8841 Human O15379
Synonyms:	HD3; RPD3; RPD3-2; SMAP45



[View online »](#)