

Product datasheet for **TA398607**

TMEM16A (ANO1) Rabbit Monoclonal Antibody [Clone ID: RM440]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM440
Applications:	IHC
Recommended Dilution:	IHC: 1:100 – 1:400 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant human DOG1 protein
Specificity:	This antibody reacts to ANO1 / DOG1 / TMEM16A
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	anoctamin 1
Database Link:	Entrez Gene 55107 Human Q5XXA6
Synonyms:	DOG1; FLJ10261; ORAOV2; TAOS2; TMEM16A



[View online »](#)