

Product datasheet for **TA398609**

Caspase 8 (CASP8) Rabbit Monoclonal Antibody [Clone ID: RM442]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM442
Applications:	WB
Recommended Dilution:	WB: 1:1000 – 1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of caspase 8
Specificity:	This antibody reacts to cleaved caspase 8
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	caspase 8
Database Link:	Entrez Gene 841 Human Q14790
Synonyms:	ALPS2B; CAP4; Casp-8; FLICE; FLJ17672; MACH; MCH5; MGC78473; OTTHUMP00000163720



[View online »](#)