

## Product datasheet for **TA398611**

### alpha Tubulin (TUBA1A) Rabbit Monoclonal Antibody [Clone ID: RM447]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM447
Applications:	ICC, WB
Recommended Dilution:	WB: 1:1000 -1:2000 dilution ICC/IF: 1:100 - 1:200 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to dC2-alpha-Tubulin
Specificity:	This antibody specifically reacts to dC2-alpha-Tubulin, a C-terminally trimmed version of alpha-tubulin lacks a carboxy-terminal glutamyl-tyrosine group.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	tubulin alpha 1a
Database Link:	<a href="#">Entrez Gene 7846 Human Q71U36</a>
Synonyms:	B-ALPHA-1; FLJ25113; hum-a-tub1; hum-a-tub2; LIS3; TUBA3



[View online »](#)