

Product datasheet for **TA398613**

PDPN Rabbit Monoclonal Antibody [Clone ID: RM449]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM449
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 -1:200 dilution WB: 1:100 - 1:500 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant human Podoplanin protein
Specificity:	This antibody reacts to human Podoplanin (Gp36)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	podoplanin
Database Link:	Entrez Gene 10630 Human Q86YL7
Synonyms:	aggrus; GP36; Gp38; GP40; HT1A-1; hT1alpha-1; hT1alpha-2; OTS8; PA2.26; podoplanin; T1-alpha; T1A; T1A-2



[View online »](#)