

## Product datasheet for **TA398614**

### **Fibrillin 1 (FBN1) Rabbit Monoclonal Antibody [Clone ID: RM463]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	RM463
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 – 1:200 dilution WB: 1:1000 - 1:25000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant full-length human Asprosin protein
Specificity:	This antibody reacts to human Asprosin.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	fibrillin 1
Database Link:	<a href="#">Entrez Gene 2200 Human P35555</a>
Synonyms:	FBN; MASS; MFS1; OCTD; SGS; WMS



[View online »](#)