

Product datasheet for **TA398615**

STAT3 Rabbit Monoclonal Antibody [Clone ID: RM465]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM465
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:50 -1:200 dilution WB: 1:1000 - 1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human Stat3
Specificity:	This antibody reacts to human Stat3. This antibody may also react to mouse or rat Stat3, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	signal transducer and activator of transcription 3
Database Link:	Entrez Gene 6774 Human P40763
Synonyms:	APRF; FLJ20882; HIES; MGC16063



[View online »](#)