

Product datasheet for **TA398619**

FAK (PTK2) Rabbit Monoclonal Antibody [Clone ID: RM472]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM472
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 -1:200 dilution WB: 1:100 - 1:400 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to residues surrounding Tyr397 of human FAK (Focal adhesion kinase 1)
Specificity:	This antibody reacts to human FAK (Focal adhesion kinase 1). This antibody may also react to mouse or rat FAK, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	protein tyrosine kinase 2
Database Link:	Entrez Gene 5747 Human Q05397
Synonyms:	FADK; FAK; FAK1; FRNK; pp125FAK



[View online »](#)