

## Product datasheet for **TA398620**

### STAT6 Rabbit Monoclonal Antibody [Clone ID: RM473]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM473
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 -1:200 dilution WB: 1:1000 - 1:2500 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to residues near the C-terminus of human Stat6
Specificity:	This antibody reacts to human Stat6.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	signal transducer and activator of transcription 6
Database Link:	<a href="#">Entrez Gene 6778 Human P42226</a>
Synonyms:	D12S1644; IL-4-STAT; STAT6B; STAT6C



[View online »](#)