

## Product datasheet for **TA398621**

### MSH2 Rabbit Monoclonal Antibody [Clone ID: RM478]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM478
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 -1:200 dilution WB: 1:1000 - 1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein within 500 amino acids of the C-terminus of human MSH2
Specificity:	This antibody reacts to human Msh2.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	mutS homolog 2
Database Link:	<a href="#">Entrez Gene 4436 Human P43246</a>
Synonyms:	COCA1; FCC1; HNPCC; HNPCC1; LCFS2



[View online »](#)