

Product datasheet for **TA398628**

EMA (MUC1) Rabbit Monoclonal Antibody [Clone ID: RM489]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM489
Applications:	WB
Recommended Dilution:	WB: 1:1000 - 1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the cytoplasmic domain of human Mucin-1
Specificity:	This antibody reacts to human Muc-1 (Mucin-1).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	mucin 1, cell surface associated
Database Link:	Entrez Gene 4582 Human P15941
Synonyms:	CD227; EMA; episialin; H23AG; KL-6; MAM6; MUC-1; MUC-1/SEC; MUC-1/X; MUC1/ZD; PEM; PEMT; PUM



[View online »](#)