

Product datasheet for **TA398629**

TIA1 Rabbit Monoclonal Antibody [Clone ID: RM490]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM490
Applications:	WB
Recommended Dilution:	WB: 1:1000 - 1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human TIA1
Specificity:	This antibody reacts to human TIA1 (Cytotoxic granule associated RNA binding protein TIA1) (T-cell-restricted intracellular antigen-1).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	TIA1 cytotoxic granule-associated RNA binding protein
Database Link:	Entrez Gene 7072 Human P31483
Synonyms:	p40-TIA-1



[View online »](#)