

## Product datasheet for **TA398630**

### alpha Synuclein (SNCA) Rabbit Monoclonal Antibody [Clone ID: RM492]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM492
Applications:	WB
Recommended Dilution:	WB: 1:1000 - 1:2000 dilution
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to residues near the C-terminus of human Alpha-synuclein
Specificity:	This antibody reacts to human, mouse or rat Alpha-synuclein (SYUA)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	synuclein alpha
Database Link:	<a href="#">Entrez Gene 6622 Human P37840</a>
Synonyms:	alpha-synuclein; MGC110988; NACP; OTTHUMP00000161559; OTTHUMP00000161561; PARK1; PARK4; PD1



[View online »](#)