

Product datasheet for **TA398633**

CDX2 Rabbit Monoclonal Antibody [Clone ID: RM495]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM495
Applications:	IHC
Recommended Dilution:	IHC: 1:100 -1:200 dilutio
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant human CDX2 protein
Specificity:	This antibody reacts to human CDX2
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	caudal type homeobox 2
Database Link:	Entrez Gene 1045 Human Q99626
Synonyms:	CDX-3; CDX3



[View online »](#)