

## Product datasheet for **TA398636**

### **LAIR1 Rabbit Monoclonal Antibody [Clone ID: RM498]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	RM498
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 -1:500 dilution WB: 1:1000 - 1:2500 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human LAIR-1 extracellular domain.
Specificity:	This antibody reacts to human LAIR-1
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	leukocyte associated immunoglobulin like receptor 1
Database Link:	<a href="#">Entrez Gene 3903 Human Q6GTX8</a>
Synonyms:	CD305; hLAIR1; LAIR-1



[View online »](#)