

## Product datasheet for **TA398637**

### galectin 9 (LGALS9) Rabbit Monoclonal Antibody [Clone ID: RM499]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM499
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 -1:200 dilution WB: 1:250 - 1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human Gal-9
Specificity:	This antibody reacts to human Gal-9 (Galectin-9)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	galectin 9
Database Link:	<a href="#">Entrez Gene 3965 Human</a> <a href="#">O00182</a>
Synonyms:	ecalectin; galectin-9; HOM-HD-21; HUAT; LGALS9A; MGC117375; MGC125973; MGC125974



[View online »](#)