

Product datasheet for TA398643

MERTK Rabbit Monoclonal Antibody [Clone ID: RM506]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM506
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 -1:200 dilution WB: 1:1000 - 1:2000 dilution
Reactivity:	Human
Host:	Rabbit
lsotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human MerTK extracellular domain.
Specificity:	This antibody reacts to human MerTK (Tyrosine-protein kinase Mer).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	MER proto-oncogene, tyrosine kinase
Database Link:	<u>Entrez Gene 10461 Human</u> <u>Q12866</u>
Synonyms:	c-mer; MER; MGC133349; RP38



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn