

Product datasheet for **TA398644**

Uroplakin II (UPK2) Rabbit Monoclonal Antibody [Clone ID: RM507]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM507
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 -1:200 dilution WB: 1:100 - 1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to residues near the C-terminus of human UPII (Uroplakin-2)
Specificity:	This antibody reacts to human UPII (Uroplakin-2)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	uroplakin 2
Database Link:	Entrez Gene 7379 Human O00526
Synonyms:	MGC138598; UP2; UPII; uroplakin-2



[View online »](#)