

Product datasheet for **TA398645**

Rabbit Monoclonal Antibody [Clone ID: RM420]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM420
Recommended Dilution:	Use at an assay dependent concentration. Assay dependen
Reactivity:	SARS-CoV
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A recombinant 2019-nCov Nucleocapsid Protein expressed in E.coli.
Specificity:	This antibody reacts to Nucleocapsid Protein of SARS-CoV-2 (2019-nCov)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)