

## Product datasheet for **TA398649**

### Rabbit Monoclonal Antibody [Clone ID: RM412]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM412
Applications:	ELISA
Recommended Dilution:	ELISA: 0.05 ug/mL – 5 ug/m
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Proteins purified from Helicobacter pylori
Specificity:	This antibody reacts to urease enzyme of Helicobacter pylori of Helicobacter pylori
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1.0 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Database Link:	-



[View online »](#)