

Product datasheet for **TA398653**

Rabbit Monoclonal Antibody [Clone ID: RM428]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM428
Recommended Dilution:	ELISA: 0.5ug/mL – 2ug/m
Reactivity:	Ape
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	25-OH Vitamin D3 Antibody-antigen Complex
Specificity:	This antibody reacts to 25-OH Vitamin D3/RM3 or 25-OH Vitamin D3/ RMH04 Complexes. It does not cross-react with Vitamin D3, Vitamin D2, and 25-OH Vitamin D2. It also does not recognize free 25-OH Vitamin D3, RM3 or RMH04.
Formulation:	50% Glycerol/PBS and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)