

Product datasheet for **TA398661**

Ighg3 Rabbit Monoclonal Antibody [Clone ID: RM218]

Product data:

| | |
|-----------------------|--|
| Product Type: | Primary Antibodies |
| Clone Name: | RM218 |
| Applications: | ELISA, WB |
| Recommended Dilution: | ELISA: 0.005 ug/mL-0.2 ug/mL WB (WB): 0.1 ug/mL-0.5 ug/mL |
| Reactivity: | Mouse |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Monoclonal |
| Immunogen: | Mouse IgG3 |
| Specificity: | This antibody reacts to the Fab region of mouse IgG3. No cross reactivity with mouse IgG1, IgG2a, IgG2b, IgM, IgA, IgE, human IgG, or rat IgG. |
| Formulation: | 50% Glycerol/PBS with 1% BSA and 0.09% sodium azide |
| Concentration: | 1 mg/mL - lot specific |
| Purification: | Protein A affinity purified from an animal origin-free culture supernatant |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |
| Gene Name: | Immunoglobulin heavy constant gamma 3 |
| Database Link: | Entrez Gene 380795 Mouse P03987 |

[View online »](#)