

## Product datasheet for **TA398663**

### Ighm Rabbit Monoclonal Antibody [Clone ID: RM109]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM109
Applications:	ELISA
Recommended Dilution:	ELISA: 0.005 ug/mL – 0.2 ug/mL
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Mouse IgM
Specificity:	This antibody reacts to mouse IgM. No cross reactivity with mouse IgG1, IgG2a, IgG2b, IgG3, IgA, IgE, human IgG, rat IgG, or goat IgG.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	immunoglobulin heavy constant mu
Database Link:	<a href="#">Entrez Gene 16019 Mouse P01872</a>



[View online »](#)