

Product datasheet for **TA398666**

Goat Monoclonal Antibody [Clone ID: RMG02]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RMG02
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 0.01 ug/mL – 0.5 ug/mL WB (WB) (nonreduced): 0.1 ug/mL-0.5 ug/mL
Reactivity:	Rabbit
Host:	Goat
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Rabbit IgG
Specificity:	This antibody reacts to the Fc region of Rabbit IgG. No cross reactivity with human IgG, rat IgG, or mouse IgG.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein G affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)